



REV	DESCRIPTION	DATE
X	X	00/00/2018



M6 Tapped holes on 25mm grid



4, 6.5mm Hole @100x100mm center

4, 6.5mm Hole @100x100mm center

Minimum Height : 150
Maximum Height : 200

GEN- Tolerances				Drawn	Material	Quantity	Finish	Date	Scale	Page No.
Linear		Angular		SKS	Aluminum-6061	01	Anodised	30/10/2020	NA	01
X	±0.5	X°	±0.5°	Checked SKS	HOLMARC OPTO-MECHATRONICS (P) Ltd	Title Motorized Vertical stage				
X.X	±0.1	All dims are in mm								
X.XX	±0.02				Approved	Part Number	Part Name	Drawing/ Model Number	Part Revision	
X.XXX	±0.005				Emc	0	Assembly	MVTS-125-50	0	