


Model No:	Carriage Size	Apreture Size	Travel	H	L	W	L1	W1
MXYTSW175175	175 x 175	75 x 75	50,50	65	175	175	75	75
MXYTSW225225	225 x 225	125 x 125	100 100	65	225	225	125	125
MXYTSW325225	325 x 225	225 x 125	200 100	75	325	225	225	125
MXYTSW325295	325 x 295	225 x 195	200 170	75	325	225	195	175
MXYTSW425295	425 x 295	325 x 195	300 170	75	425	225	195	175

"INFORMATION ONLY NOT FOR MANUFACTURING"

ALL DIMENSIONS ARE IN MM
UNLESS OTHERWISE STATED



HOLMARC
OPTO-MECHATRONICS PVT. LTD

MXYTSW Series XY Stage

THIS INFORMATION IS CONFIDENTIAL AND PROPRIETARY TO HOLMARC, ITS WORLDWIDE SUBSIDIARIES AND/OR AFFILIATES. IT MAY NOT BE DISCLOSED TO ANYONE, OTHER THAN HOLMARC PERSONNEL, WITHOUT WRITTEN AUTHORIZATION FROM AN AUTHORIZED REPRESENTATIVE OF HOLMARC.