



Model No.	Carriage Size	Travel	Height	L	W
MXYTS150150	150x150	50,50	65	150	150
MXYTS170170	170x170	100,100	65	170	170
MXYTS270170	270x170	200,100	75	270	170
MXYTS270240	270x240	200,170	75	270	240
MXYTS370240	370x240	300,170	75	370	240

"INFORMATION ONLY NOT FOR MANUFACTURING"

ALL DIMENSIONS ARE IN MM
UNLESS OTHERWISE STATED

Scale : NA



XY MOVING PLATFORM

THIS INFORMATION IS CONFIDENTIAL AND PROPRIETARY TO HOLMARC, ITS WORLDWIDE SUBSIDIARIES AND/OR AFFILIATES. IT MAY NOT BE DISCLOSED TO ANYONE, OTHER THAN HOLMARC PERSONNEL, WITHOUT WRITTEN AUTHORIZATION FROM AN AUTHORIZED REPRESENTATIVE OF HOLMARC