



Plunger knob



Model No:	Height-H(mm)
PH-16D-30	30
PH-16D-40	40
PH-16D-50	50
PH-16D-75	75
PH-16D-100	100
PH-16D-150	150

"INFORMATION ONLY NOT FOR MANUFACTURING"

ALL DIMENSIONS ARE IN MM  
UNLESS OTHERWISE STATED



**HOLMARC**  
OPTO-MECHATRONICS PVT. LTD

PH-16D Series Post  
Holder

THIS INFORMATION IS CONFIDENTIAL AND PROPRIETARY TO HOLMARC , ITS WORLDWIDE SUBSIDIARIES AND/OR AFFILIATES.IT MAY NOT BE DISCLOSED TO ANYONE, OTHER THAN HOLMARC PERSONNEL, WITHOUT WRITTEN AUTHORIZATION FROM AN AUTHORIZED REPRESENTATIVE OF HOLMARC