


Model No:	Type	Suitable Post Series	Size (mm)
PB - 6550	Square-Type	12 & 16	65 x 50 x 9
MPB - 6550	Square-Type	12 & 16	65 x 50 x 9

"INFORMATION ONLY NOT FOR MANUFACTURING"

ALL DIMENSIONS ARE IN MM  
UNLESS OTHERWISE STATED



**HOLMARC**  
OPTO-MECHATRONICS PVT. LTD

Post Bases

THIS INFORMATION IS CONFIDENTIAL AND PROPRIETARY TO HOLMARC , ITS WORLDWIDE SUBSIDIARIES AND/OR AFFILIATES.IT MAY NOT BE DISCLOSED TO ANYONE, OTHER THAN HOLMARC PERSONNEL, WITHOUT WRITTEN AUTHORIZATION FROM AN AUTHORIZED REPRESENTATIVE OF HOLMARC