



"INFORMATION ONLY NOT FOR MANUFACTURING"

Model No:	Carriage Size	Aperture Size	Travel	H	L	W	L1	W1
XYTSW175175	175 x 175	75 x 75	50,50	65	175	175	75	75
XYTSW225225	225 x 225	125 x 125	100100	65	225	225	125	125
XYTSW325325	325 x 225	225 x 125	200100	75	325	225	225	125
XYTSW325295	325 x 295	225 x 195	200170	75	325	225	195	175
XYTSW425295	425 x 295	325 x 195	300170	75	425	425	195	175

ALL DIMENSIONS ARE IN MM
UNLESS OTHERWISE STATED



XYTSW Series

THIS INFORMATION IS CONFIDENTIAL AND PROPRIETARY TO HOLMARC, ITS WORLDWIDE SUBSIDIARIES AND/OR AFFILIATES. IT MAY NOT BE DISCLOSED TO ANYONE, OTHER THAN HOLMARC PERSONNEL, WITHOUT WRITTEN AUTHORIZATION FROM AN AUTHORIZED REPRESENTATIVE OF HOLMARC