



Model No:	Carriage Size	Travel	Height (mm)	L(mm)	B(mm)
XYTS150150	150x150	50,50	65	150	150
XYTS170170	170x170	100100	65	170	170
XYTS270170	270x170	200100	75	270	170
XYTS270240	270x240	200170	75	270	240
XYTS370240	370x240	300170	75	370	240

"INFORMATION ONLY NOT FOR MANUFACTURING"

ALL DIMENSIONS ARE IN MM
UNLESS OTHERWISE STATED

Scale : NA



XYTS Series

THIS INFORMATION IS CONFIDENTIAL AND PROPRIETARY TO HOLMARC, ITS WORLDWIDE SUBSIDIARIES AND/OR AFFILIATES. IT MAY NOT BE DISCLOSED TO ANYONE, OTHER THAN HOLMARC PERSONNEL, WITHOUT WRITTEN AUTHORIZATION FROM AN AUTHORIZED REPRESENTATIVE OF HOLMARC.