



"INFORMATION ONLY NOT FOR MANUFACTURING"

ALL DIMENSIONS ARE IN MM  
UNLESS OTHERWISE STATED

Scale : NA



**HOLMARC**  
OPTO-MECHATRONICS PVT. LTD

Model No:	Drive	Travel
TC XZ-M $\mu$ 10	Micrometer 10 $\mu$	5mm
TS XZ-D1FM1	Differential Micrometer 1 $\mu$	2.5mm

XZ Positioner

THIS INFORMATION IS CONFIDENTIAL AND PROPRIETARY TO HOLMARC, ITS WORLDWIDE SUBSIDIARIES AND/OR AFFILIATES. IT MAY NOT BE DISCLOSED TO ANYONE, OTHER THAN HOLMARC PERSONNEL, WITHOUT WRITTEN AUTHORIZATION FROM AN AUTHORIZED REPRESENTATIVE OF HOLMARC