



Model No:	Maximum Size	A	H	T
ALM - B - 50	50	85	90	15
ALM B - 90	90	126	145	18
ALM B - 135	135	165	200	25
ALM B -165	165	200	245	25

"INFORMATION ONLY NOT FOR MANUFACTURING"

ALL DIMENSIONS ARE IN MM  
UNLESS OTHERWISE STATED



**HOLMARC**  
OPTO-MECHATRONICS PVT. LTD

Bar Type Adjustable  
Lens Mount

THIS INFORMATION IS CONFIDENTIAL AND PROPRIETARY TO HOLMARC , ITS WORLDWIDE SUBSIDIARIES AND/OR AFFILIATES.IT MAY NOT BE DISCLOSED TO ANYONE, OTHER THAN HOLMARC PERSONNEL, WITHOUT WRITTEN AUTHORIZATION FROM AN AUTHORIZED REPRESENTATIVE OF HOLMARC