


Model No:	Optic Diameter : D	Filter Cell Model	D1	D2	D3
QC - FCH - 12.5	15	34	38	18	10
QC - FCH - 25	25	44	50	30	22
QC - FCH - 38	30	54	62	44	34
QC - FCH - 50	50	79	70	54	42

"INFORMATION ONLY NOT FOR MANUFACTURING"

ALL DIMENSIONS ARE IN MM
UNLESS OTHERWISE STATED



HOLMARC
OPTO-MECHATRONICS PVT. LTD

Quick Change Multiple
Filter Holder

THIS INFORMATION IS CONFIDENTIAL AND PROPRIETARY TO HOLMARC , ITS WORLDWIDE SUBSIDIARIES AND/OR AFFILIATES.IT MAY NOT BE DISCLOSED TO ANYONE, OTHER THAN HOLMARC PERSONNEL, WITHOUT WRITTEN AUTHORIZATION FROM AN AUTHORIZED REPRESENTATIVE OF HOLMARC