


Model No:	Optic Size	D	D1	L	T
XYZLM-25	25/1"	25	25	65	20
XYZLM-50	50/2"	50	50	90	20

"INFORMATION ONLY NOT FOR MANUFACTURING"

ALL DIMENSIONS ARE IN MM  
UNLESS OTHERWISE STATED



**HOLMARC**  
OPTO-MECHATRONICS PVT. LTD

XYZ Lens Mounts

THIS INFORMATION IS CONFIDENTIAL AND PROPRIETARY TO HOLMARC, ITS WORLDWIDE SUBSIDIARIES AND/OR AFFILIATES. IT MAY NOT BE DISCLOSED TO ANYONE, OTHER THAN HOLMARC PERSONNEL, WITHOUT WRITTEN AUTHORIZATION FROM AN AUTHORIZED REPRESENTATIVE OF HOLMARC