



CIN : U33125KL1993PTC006984
GST No: 32AAACH9492C1ZQ

B-7., H.M.T. Industrial Estate, H.M.T. PO, Kalamassery
Kochi - 683 503 Kerala, India.
Email : sales@holmarc.com, PH : 91-484-2540075



QA/Test Report

Equipment name : Silar coating unit

Date:


Model No./ Serial no. :

Customer/Job card No.

Standard Test Procedure: (

All the tests are performed based on the theory of the respective experiment and/or by calibrated instruments.)

SI No	Specific Test Conducted	Test Result	Specification/ Accepted values	Pass /Fail	Remarks
1.	Basic Inspection		Model No., Painting , Finishing , stability and other visual parameters		
2.	Vibration and noise checking		Run a sample program and check		
3.	Keyboard/ Touch screen checking		Enter all the alpha numeric values and check its functioning		
4.	Up-Down travel and limit checking		Move the sample holder up side and down side till it reaches the limit		
5.	Travel length		Verify the travel length by running a sample program with specific travel length say 50 mm		
6.	Travel speed		Verify the travel speed by running a sample program with specific travel speed (both Dip and retrieval) in um/Sec and mm/Sec		
7.	Dip & dry duration		Run a sample program with specific Dip and Dry duration. Sec., Min. & Hour		
8.	Sample holder positioning. (centering to the beaker)		Run a sample program with Dip in all the beakers and check the centering at each dip position.		
9.	Magnetic stirrer Function		Check the magnetic stirrer with beaker and run @ max rpm.		
10.	Heater Function		Turn on the heater and check the temperature.		
11	PC -Software testing		Compose and run a sample program in PC-Software.		



Checked by

Approved by

JOLLY CYRIAC
Managing Director



Mailing Address
B-7., H.M.T. Industrial Estate,
H.M.T. PO, Kalamassery
Kochi - 683 503 Kerala, India.



E-mail Contact
sales@holmarc.com



Phone Contact
Tel: 91-484-2540075



www.holmarc.com

Business Hours
Monday to Saturday: 9:00am - 5:00pm

