



CIN : U33125KL1993PTC006984
GST No: 32AAACH9492C1ZQ

ISO 9001:2015

HOLMARC
OPTO-MECHATRONICS PVT.LTD

B-7., H.M.T. Industrial Estate, H.M.T. PO, Kalamassery
Kochi - 683 503 Kerala, India.
Email : sales@holmarc.com, PH : 91-484-2540075



QA/Test Report

Equipment name : Syringe pump.

Date:


Model No./ Serial no. :

Customer/Job card No.

Standard Test Procedure: (

All the tests are performed based on the theory of the respective experiment and/or by calibrated instruments.)

Sl No	Specific Test Conducted	Test Result	Specification/ Accepted values					Pass/Fail
1.	Basic Inspection		Model No., Painting , Finishing , stability and other visual parameters					
2.	Vibration and noise checking		Run a sample program and check					
3.	Keyboard/ Touch screen checking		Enter all the alpha numeric values and check its functioning					
4.	Limit switch positioning		Adjust the far limit sensor position as suitable for 20ml glass syringe.					
5.	Limit checking		Run home and far -check the stage movement stops in home and far position.					
6.	Syringe pump		Run sample programs and check the flow rates and corresponding volumes.					
7.	PC -Software testing		Run sample programs in software also					
	Flow rate		Sl. No.	Set(ml/min)	Measured(ml/min)	Error (%)	Permissible Error (%)	
			1					
			2					
			3					
	Volume		1					
			2					
			3					



[Signature]

Checked by

.....

Approved by

.....

JOLLY CYRIAC
Managing Director



Mailing Address
B-7., H.M.T. Industrial Estate,
H.M.T. PO, Kalamassery
Kochi - 683 503 Kerala, India.



E-mail Contact
sales@holmarc.com



Phone Contact
Tel: 91-484-2540075



www.holmarc.com

Business Hours
Monday to Saturday: 9:00am - 5:00pm

