



CIN : U33125KL1993PTC006984  
GST No: 32AAACH9492C1ZQ

ISO 9001:2015

**HOLMARC**  
OPTO-MECHATRONICS PVT.LTD

B-7., H.M.T. Industrial Estate, H.M.T. PO, Kalamassery  
Kochi - 683 503 Kerala, India.  
Email : sales@holmarc.com, PH : 91-484-2540075



### QA/Test Report

Equipment name : Honeycomb Breadboard Date:

Model No./ Serial no. :

Customer/Job card No. ....

Standard Test Procedure: (

All the tests are performed based on the theory of the respective experiment and/or by calibrated instruments.)


SI No	Specific Test Conducted	Test Result	Specification/ Accepted values	Pass/Fail
1.	Basic Inspection		Model No., Painting , Finishing , stability and other visual parameters	
2.	Vibration and noise checking		Assemble the bread board and study the vibration damping	
3.	Dimension		Check all the dimensions using the scale	
4.	Compressor		Switch on the compressor and check the pressure	
5.	Air leakage testing		Pressurize the supports and keep for long time and check the pressure changed or not	

Checked by

.....

Approved by

.....



**JOLLY CYRIAC**  
Managing Director



Mailing Address  
B-7., H.M.T. Industrial Estate,  
H.M.T. PO, Kalamassery  
Kochi - 683 503 Kerala, India.



E-mail Contact  
sales@holmarc.com



Phone Contact  
Tel: 91-484-2540075



www.holmarc.com

Business Hours  
Monday to Saturday: 9:00am - 5:00pm

Connect with us on

