



**HOLMARC**  
OPTO-MECHATRONICS PVT.LTD



### Quality Test Plan

**Equipment name** : Bread board

**Mode/ Serial no** : BB series

Sl No	Specific Test Conducted	Test method	No. of times repeated
1.	Basic Inspection	Model No., Painting , Finishing , stability and other visual parameters	1
2.	Vibration and noise checking	Assemble the bread board and study the vibration damping	1
3.	Dimension	Check all the dimensions using the scale	3
4.	Compressor	Switch on the compressor and check the pressure	1
5.	Air leakage testing	Pressurize the supports and keep for long time and check the pressure changed or not	3