

भारत सरकार
परमाणु ऊर्जा विभाग
क्रय एवं भण्डार निदेशालय
इंदौर क्षेत्रीय क्रय एवं भण्डार इकाई
राजा रामन्ना प्रगत प्रौद्योगिकी केन्द्र



फार्म संख्या डीपीएस-पी-23 / Form No.: DPS P-23
Government of India
Department of Atomic Energy
Directorate of Purchase & Stores
Indore Regional Purchase & Stores Unit

Raja Ramanna Centre for Advanced Technology
PO : CAT, INDORE - 452013 (M.P.) INDIA

पो आ : केट. इंदौर 452013 (म प) भारत
दूरभाष / Tel. #: 091-731-248-8795/8786/8797/8798, फेक्स / Fax#: 091-731- 2321345, 2488800, वेबसाइट / Website : www.rrcat.gov.in

RRCAT GSTIN : 23AAAGR0577G1ZF

क्रय आदेश सं

क्रय आदेश / PURCHASE ORDER

Date: 08 FEB 2019

PURCHASE ORDER NO.: DPS/IRPU/PG-II/61024/MES/PO/57550

Your Offer No.:

Ref.:HOMPL/R&D/QTN/2018-19/
793-A dtd 12.09.2018

To,
HOLMARC OPTO-MECHATRONICS PVT LTD
B-7 H.M.T. INDUSTRIAL ESTATE,
HMT PO., KALAMASSERY
KOCHI (KERALA)-683503 PIN : 683503
PHONE : 0484-2540075 FAX : 0484-2540882
E-mail : sales@holmarc.com

08 NOV 2019

प्रिय महोदय,
निम्नलिखित सामान की पूर्ति के लिए दिनांक 14-AUG-18 की निविदा सं DPS/IRPU/PG-II/61024 में दिया गया आपका निविदा भारत के राष्ट्रपति की ओर से, उक्त निविदा (फार्म संख्या डीपीएस-पी-101/102 जैसा लागू हो) के पिछले पृष्ठ पर दी गई शर्तों के अनुसार तथा नीचे विनिर्दिष्ट सामान की मात्रा तथा उसकी कीमत की सीमा तक इस क्रय आदेश में दी गयी शर्तों के अनुसार, स्वीकार किया जाता है।

Dear Sirs,

Your offer contained in the Tender No. DPS/IRPU/PG-II/61024 dated 14-AUG-18 to supply the under mentioned stores is accepted on behalf of the President of India, subject to the terms & conditions of contract printed at the back of the said tender (Form No. D.P.S.-P.-101/102 as applicable) and subject to the terms and conditions contained in this Purchase Order to the extent of the quantity and at the price specified below :

क्रम सं S.No.	विवरण / Description	मात्रा Quantity	इकाई Unit	दर इकाई Unit Rate (Rs.)	कुल मूल्य Total Value (Rs.)
1.	Supply, Installation, Commissioning & Demonstration of Spin Coater Compatible with Glove Box Model : HO-TH-05G Make : Holmarc Optomechatronics (Detailed specifications as per Annexure-I encl.)	2 Nos.	Each	194450.00	388900.00
Grand Total of all the Items:					388900.00

Other Charges:

- + Packing & Forwarding Charges extra at actual but not to exceed @ 2% of basic price.
- + Goods & Service Tax (GST) where applicable will be paid extra at actual during the delivery period stipulated in the Purchaser Order. In terms of notifications issued by the Central Government & State Government, R&D units of Department of Atomic Energy are entitled for IGST @ 5% as applicable for stores covered under the Purchase Order.
- It would be the responsibility of the contractor to ensure that relevant certificate is obtained from the purchaser before affecting the delivery of goods ordered failing which the excess tax paid by the contractor shall not be reimbursed by the purchaser. The contractor should furnish a certificate along with Invoice as per Appendix 'A' enclosed.
- + Freight charges extra at actuals on production of documentary evidence basis but it should not exceed @ 1% of basic price.
- + T.D.S. as applicable shall be deducted at source on payment made to the supplier of taxable goods.

Notes:

- Installation & Commissioning of the supplied material shall be carried out by the supplier at our RRCAT Site, Indore without any extra cost.

मूल्य / Prices	Free Delivery at Stores, DPS, RRCAT, INDORE-13	मुपुर्दगी अवधि / Delivery Period	on or before 08 JAN 2019
मुपुर्दगी का स्थान / Place of Delivery	Asst. Stores Officer, DPS, RRCAT, Indore-13		

कृपया यह नोट करें कि उपरिनिर्दिष्ट निविदा फार्म संख्या डीपीएस-पी-101/102 (जैसा लागू हो) के पिछले पृष्ठ पर तथा इस कार्य आदेश में दी गई शर्तों ही ठेके की एकमात्र आधार होगी।

Please note that the terms and conditions printed at the back of the tender form No. D.P.S.-P.-101/102 (as applicable) referred to above and this Purchase Order shall be the sole repository of the contract.

कृपया यहाँ मंगलन फार्म में तुरंत प्राप्ति सूचना दें।

Please acknowledge receipt immediately in the form sent herewith.

भवदीय / Your faithfully,

(G. RAMU)

ASSTT. PURCHASE OFFICER

भारत के राष्ट्रपति के लिए और उनकी ओर से (केता)

For and on behalf of the President of India (The Purchaser)

SUPPLIER COPY

Annexure-I

Purchase Order No. : DPS/IRPU/PG-II/61024/MES/PO/57550

● Technical Specifications:

Item name: Spin coater compatible with glove box

Quantity: 02 Nos.

Model: HO-TH-05G

1. Design : Designed to use inside a glove box, A separate controller with user friendly front panel having keyboard and LCD display to operate the spin coater from remote location.
2. Spin chamber :
 - i. Material: Nylon, flat lid with reduced height easily installable inside a glove box
 - ii. Foot print area: ~ 300 mm x 300 mm
 - iii. Top: transparent safety lid to protect the surroundings by keeping the spin off solution inside.
 - iv. The entire Spin chamber unit will be placed inside the glove box without fixing it on the base of the glove box, Therefore it should have vibration resistant design for smooth operation.
3. Actuator : Precision brush less DC servo motor with accurate speed and acceleration control
4. Spinning speed : ~ 60- 9999 RPM